

Respiratory Panel Real Time PCR Detection

SUMMARY AND EXPLANATION OF THE TEST

The respiratory panel real time PCR detection reagent is designed for simultaneous detection of the following pathogens

Strep A	Bordetella pertussis
RNase P	Chlamydia pneumonia
Flu A	Mycoplasma pneumonia
Flu B	Adenovirus
RSV	Legionella pneumophila
Enterovirus	Klebsiella pneumonia group
HMPV	Staphylococcus aureus
Rhinovirus	Moraxella catarrhalis
PIV4	Bordetella parapert
PIV3	MRSA
PIV2	H Bocavirus
PIV1	Haemophilus influenzae
229E	SC2-N1
HKU1	SC2-N2
NL63	Haemophilus influenzae B
OC43	Streptococcus pneumoniae

REAGENTS AND MATERIALS PROVIDED

Item	Component	Quantity		Volume /T
		100 Tests / Box	500 Tests / Box	
1	PCR Mix 1-8	1.5 mL/ tube each	1.5 mL/ tube × 5 each	14 µL
2	Enzyme Mix	850 µL/ tube	1550 µL/ tube	1 µL
3	Positive Control	100 µL/ tube × 1	250 µL/ tube	5 µL
4	Negative Control	100 µL/ tube × 1	250 µL/ tube	5 µL

MATERIALS MAY REQUIRED BUT NOT PROVIDED

1. Real-Time PCR Thermocycler compatible with JOE ([VIC or HEX is equivalent to JOE](#)), ROX, FAM and Cy5 channels.
2. Specimen Collection Containers
3. RNA extraction reagent
4. PCR tubes compatible with Thermocycler.
5. 1.5 mL centrifuge tubes (RNase-Free)
6. Vortex mixer
7. High-speed centrifuge
8. Micropipette (0.5-2 µL, 1-10 µL)
9. Pipette tips with filters
10. Personal protective material

WARNING AND PRECAUTIONS

1. For research use only. Not for diagnostic procedures.
2. Check the expiration date.
3. Follow Good Laboratory Practices: wear appropriate protective clothing and use disposal gloves and protective eyewear. Do not eat, drink, or smoke in designated work areas. Wash hands thoroughly after handling specimens and kit reagents.
4. The testing workflow must be one-directional to minimize contamination risks (allocate segregated areas for each step): it should be started from RNA Purification Area, then move to the Reaction Setup Area, followed by Amplification and Detection Area. Do not bring samples, equipment, and reagents to the area in which the previous step was performed and always change gloves when changing areas.
5. Regular decontamination of commonly used equipment is recommended, especially for micropipettes and work surfaces.
6. Specimens must be treated as potentially infectious sources as well as all reagents and materials that have been exposed to the samples and handled in the same manner as an infectious agent. Take appropriate precautions during specimen collection, storage, handling, and disposal

- in accordance with the guidance from country or region authorities.
7. To ensure optimal performance of the test, always follow appropriate procedures for specimen collection, transport, storage, and processing. Improper procedures may lead to false negative results.
 8. Nucleic acids must first be extracted from specimens using a RNA Extraction Kit prior testing. Do not perform a RT-PCR Assay directly using specimens without extraction.
 9. Appropriate precautions should be taken to monitor contamination and preserve the purity of kit components and reactions. Avoid microbial and nuclease (RNase/DNase) contamination of specimens and kit components. Avoid the spread of aerosols when handling or uncapping specimens.
 10. To minimize cross-contamination, open only one tube at a time in the process.
 11. Always change pipette tips between liquid transfers. To minimize cross contamination, it's recommended to use aerosol-barrier pipette tips.
 12. Always use disposable gloves and regularly check that they are not contaminated with sample materials. Discard gloves if they become contaminated.

TRANSPORT AND STORAGE INSTRUCTIONS

Store the kit frozen at -20±5°C and avoid repeated freeze-thaw cycles (less than 5 frozen-thaw times allowed). The product is valid for 6 months, and the expiry date of the product is shown on the label.

SPECIMEN TYPE and PRESERVATION

Specimen type: oropharyngeal swabs, nasopharyngeal swabs and sputum.

Specimen preservation: process specimens for viral RNA extraction and nucleic acid detection as soon as possible. Specimens can be stored at 4°C for 24 hours, or frozen at -80°C for 1 year, if not processed immediately. Avoid repeated freeze-thaw cycles during transport and storage of the specimens.

ASSAY PCR REACTION SETUP (at reagent preparation area)

A. Reagent Preparation

To calculate the number (n) of PCR reaction tubes (n= number of specimens to be tested + number of positive controls + number of negative controls), add n x 14 µL of PCR Mix and n x 1 µL of RT-PCR Enzyme Mix to a clean tube, vortex to mix well, centrifuge and aliquot 15 µL of the mix into each PCR reaction tube.

B. Addition of Viral RNA

Add 5 µL of viral RNA extracts of the specimen or positive/negative control to each PCR tube according to the reaction preparation table below for testing. Then immediately cap the tube tightly, centrifuge the mixture and move it to the amplification area for RT-PCR detection.

PCR Reaction Preparation:

Component	Volume per test
PCR Mix	14 µL
RT-PCR Enzyme Mix	1 µL
Negative Control / Positive Control / RNA of specimen to be tested	5 µL
Total Reaction Volume	20 µL

C. RT-PCR amplification (at Amplification Area, use Applied Biosystems 7500 Real-Time PCR System as an example)

1. Place the complete PCR reaction tubes into the fluorescent quantitative PCR analyzer and label positive control, negative control and testing specimen on the Thermocycler software.
2. Enter the PCR reaction volume: 20 µL.
3. Select fluorescence detection channels: JOE, ROX, FAM and Cy5 channels.
4. Select PCR cycle parameter setting following the instructions below, save the file and start the PCR run.
5. Analyze the results according to your thermocycler manufacturer's instructions. If the positive and negative control samples do not meet their respective requirements, the PCR run must be repeated.

Note: The positive and negative control can be used directly without purification.

Cycles	Step	Time	Temperature
1	Contamination removal	2 minutes	37 °C
1	Reverse transcription	5 minutes	50 °C
1	Initial denaturation	2 minutes	95 °C
40	Denaturation	5 seconds	95 °C
	Annealing/Extension	20 seconds	60 °C

QUALITY CONTROL

Quality control requirements must be performed in conformance with local, accreditation requirements and the user's laboratory's standard quality control procedures.

channel.

DETERMINATION OF RESULTS

- The results can only be determined when the results of both controls meet their requirements.

Multiplex PCR Setup 1

Channel	Pathogens	Results Interpretation
Fam	Strep A	
JOE	RNase P	
ROX	Flu A	
Cy5	Flu B	

Multiplex PCR Setup 2

Channel	Pathogens	Results Interpretation
JOE	RSV	
Fam	Enterovirus	
ROX	HMPV	
Cy5	Rhinovirus	

Multiplex PCR Setup 3

Channel	Pathogens	Results Interpretation
JOE	PIV4	
FAM	PIV3	
ROX	PIV2	
Cy5	PIV1	

Multiplex PCR Setup 4

Channel	Pathogens	Results Interpretation
HEX	229E	
FAM	HKU1	
ROX	NL63	
Cy5	OC43	

Multiplex PCR Setup 5

Channel	Pathogens	Results Interpretation
Fam	Bordetella pertussis	
ROX	Chlamydia pneumonia	
Cy5	Mycoplasma pneumonia	
VIC	AV	

Multiplex PCR Setup 6

Channel	Bacteria	Results Interpretation
JOE	Legionella pneumophila	
Rox	<u>Klebsiella pneumonia group</u>	
FAM	Staphylococcus aureus	
Cy5	Moraxella catarrhalis	

Multiplex PCR Setup 7

Channel	Bacteria	Results Interpretation
JOE	B. paraperit	
Rox	MRSA	
FAM	H Boca	
Cy5	Haemophilus influenzae	

Multiplex PCR Setup 8

Channel	Pathogens	Results Interpretation
JOE	SC2-N1	
FAM	SC2-N2	
ROX	Haemophilus influenzae B	
Cy5	Streptococcus pneumoniae	

COMPATIBILITY OF REAL-TIME PCR THERMOCYCLERS

Manufacturer	Model
Applied Biosystems	7500 Real-Time PCR System
	7300 plus Real-Time PCR System

REFERENCE

1. National Health Commission & State Administration of Traditional Chinese Medicine. Diagnosis and Treatment Protocol for Novel Coronavirus Pneumonia (Trial Version 7). [R]. March 3, 2020

Symbols used

Symbol	Meaning
	Manufacturer
	Authorized representative in the European Community
	Catalogue number
	Batch code
	The number of tests
	Caution
	Use by date
	Consult instruction for use
	Date of manufacture
	Storage temperature limit

Version 250515



Company name: Biomiga, Inc.
 Address: 10637 Roselle Street, Suite C, San Diego, CA 92121
 USA
 Tel: 858-603-3219
 Email: info@biomiga.com



Company name: Lotus NL B.V.
 Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The
 Hague, Netherlands.
 Email: peter@lotusnl.com

Intended Use:

For research use only (RUO). Not for use in diagnostic procedures.