### **Endotoxin Removal Procedure**

This protocol is designed to remove endotoxin after the plasmid DNA is purified (Buffers can be scaled up or down accordingly).

- 1. Add 0.2 volume of EndoClean buffer to the plasmid sample in a 1.5 mL sterile tube. (For example, add 0.2 mL endoClean buffer to 1 mL plasmid sample.
- 2. Mix by vortexing the tube a few times and put on ice for about 10 minutes until the solution is clear without turbidity. (Rocking the sample in a cold room for 10 minutes is recommended if it is available).
- 3. Mix well again by inverting the tube a few times.
- **4**. Incubate the tube at 55 °C water bath for about 5 minutes and the solution shall be turbid.
- **5.** Centrifuge at top speed at **room temperature** for 5 minutes with the brake off.
- **6.** Carefully transfer the upper clear layer solution to 1.5 mL tube. Split in two 1.5 mL tubes, 400- 450 uL each tube.
- 7. Precipitate plasmid DNA with 0.1 volume of 3 M NaAc (pH 5.2) and 2 volume of 100% ethanol.

Note: The EndoClean-treated plasmid DNA is ready transfection of endotoxin sensitive cell lines.

### **Contents**

Important Notes	1
Introduction	2
Storage and Stability	2
Kit Contents	
Before Starting	3
EZgene <sup>TM</sup> Plasmid Miniprep Protocol	
Trouble Shooting Guide	
Endotoxin Removal Protocol	

## **Important Notes:**

<u>Plasmid Copy Numbers:</u> The yield of plasmid DNA depends on the origin of the replication and the size of the plasmid. The protocols are optimized for high copy number plasmid purification. For low copy number plasmids, both the culture volume and the buffer volume need to be scaled up 2 to 3 times.

Plasmid	Origin	Copy Numbers	Expected Yield (µg per 1 mL)
pSC101	pSC101	5	0.1-0.2
pACYC	P15A	10-12	0.4-0.6
pSuperCos	pMB1	10-20	0.4-1
pBR322	pMB1	15-20	0.6-1
pGEM <sup>R</sup>	Muted pMB1	300-400	6-7
pBluescript <sup>R</sup>	ColE1	300-500	6-8
pUC	Muted pMB1	500-700	8-12

Host Strains: The strains used for propagating plasmid have significant influence on yield. Host strains such as Top 10, DH5a, and C600 yield high-quality plasmid DNA. endA+ strains such as JM101, JM109, JM110, HB101, TG1 and their derivatives, normally have low plasmid yield. We recommend transform plasmid to a endA- strain if the yield is not satisfactory. For endA+ strains, we recommend use Buffer KB to remove the endonuclease before washing step.

<u>Culture Medium:</u> This procedure is designed for isolating plasmid grown in standard LB medium (Luria Bertani) to density of  $OD_{600}\,2.0$  to 3.0. If rich medium such as TB or 2xYT are used, make sure the cell density doesn't exceed 3.0 ( $OD_{600}$ ). A high ratio of cell density over lysis buffers result in low DNA yield and purity. For over amount of cell numbers, either reduce the biomass or scale up the volumes of Buffer F1, B1 and N1.

### Introduction

Key to the kit is our proprietary DNA binding systems that allow the high efficient binding of DNA to our ezBind<sup>TM</sup> matrix while proteins and other contaminates are removed under certain optimal conditions. Nucleic acids are easily eluted with sterile water or TE buffer. The purified DNA is ready for downstream applications such as cloning/subcloning, RFLP, sequencing, and transfection of HEK293 cells. The yield from 1 mL culture is typical around 8 to 12 ug. The mini column has a DNA binding capacity of 40 ug.

# **Storage and Stability**

Buffer F1 should be stored at 4°C once RNase A is added. All other materials can be stored at room temperature. The Guaranteed shelf life is 18 months from the date of purchase.

### **Kit Contents**

Catalog Number	D1218-01	D1218-02
Preps	50	250
Lysate Clearance Column	50	250
DNA Columns	50	250
Buffer F1	15 mL	60 mL
Buffer F2	12 mL	60 mL
Buffer F3	12 mL	60 mL
Buffer KB	25 mL	125 mL
DNA Wash Buffer	12 mL	50 mL
Elution Buffer	10 mL	30 mL
RNase A	50 μL	210 μL

<sup>\*</sup>Buffer F3 contains acid acid salts, wear gloves and protective eyewear when handling.

# **Trouble Shooting Guide**

Low Yield	Poor Cell lysis.	Resuspend pellet
Low Tield	1 ooi cen 1ysis.	thoroughly by
		votexing and pipeting
		prior adding Buffer
		F2.
		• Make fresh Buffer F2
		if the cap had not
		been closed tightly.
		(Buffer F2: 0.2N
		NaOH and 1%SDS).
Low Yield	Bacterial culture	Grow bacterial 12-16
	overgrown or not fresh.	hours. Spin down cultures
		and store the pellet at
		-20°C if the culture is not
		purified the same day. Do
		not store culture at 4°C
		over night.
Low Yield	Low copy-number	Increase culture volume
	plasmid.	(up to 10mL for
		Minipreps, 200mL for
		Midipreps, 400 mL for
		Maxipreps and 3L for
		Megapreps). Scale up the
		volume of buffers
		accordingly.
No DNA	Plasmid lost in Host <i>E.coli</i>	Prepare fresh culture.
Genomic DNA	Over-time incubation after	Do not vortex or mix
contamination	adding Buffer F2.	aggressively after adding
		Buffer F2. Do not incubate
		more than 5 minutes after
		adding solution B1.
RNA contamination	RNase A not added to	Add RNase A to Buffer
	Buffer F1.	F1.
Plasmid DNA floats out of	Ethanol traces not	Make sure that no ethanol
wells while running in	completely removed from	residual remaining in the
agarose gel, DNA doesn't	column.	silicon membrane before
freeze or smell of ethanol		elute the plasmid DNA.
		Re-centrifuge or vacuum
		again if necessary.
L	L	

<sup>\*\*</sup>Add 48 mL (D1218-01) or 200 mL (D1218-02) 96-100% ethanol to DNA Wash Buffer before use.

#### 11. The DNA concentration can be calculated as follows,

DNA concentration = Absorbance 260 nm x 50 x dilution factor  $(\mu g/ml)$ .

**Note:** The DNA is ready for down stream applications such as cloning, sequencing, RFLP or transfection of 293 cells.

**Note**: It's highly recommended to remove the endotoxin if the DNA is used for endotoxin-sensitive cell lines, primary cultured cells or microinjection.

Note: Use less elution buffer if high DNA concentration is desired.

# **Before Starting**

Prepare all components and get all necessary materials ready by examining this instruction booklet and become familiar with each steps.

### **Important:**

- RNase A: Spin down RNase A vial briefly. Add the RNase A solution to Buffer F1 and mix well before use.
- Add 48 mL (1211-01) or 200 mL (D1211-02) 96-100% ethanol to DNA Wash Buffer before use
- Buffer F2 precipitates below room temperature, it is critical to warm up the buffer at 37°C to dissolve the precipitates before use.
- Keep the cap tightly closed for Buffer F2 after use.
- Ensure the availability of centrifuge capable of 13,000 rpm.
- Carry out all centrifugations at room temperature.
- Buffer F3 contains acid acid, wear gloves and protective eyewear when handling.

# **Materials supplied by users:**

- 96- 100% ethanol.
- 1.5 mL microcentrifuge tubes
- High speed microcentrifuge.

# EZgene TM Plasmid Miniprep Protocol

- 1. Harvest 1-2 mL of <u>fresh</u> bacterial culture by centrifugation for 30 seconds at 12,000 rpm. Pour off the supernatant and blot the inverted tube on a paper towel to remove residue medium. Remove the residue medium completely.
- 2. Add 200 μL Buffer F1 and completely resuspend bacterial pellet by vortexing or pipetting (Complete resuspension is critical for optimal yields).
- 3. Add 200 µL Buffer F2, mix by inverting 10 times (**do not vortex**) and incubate at room temperature for 2 minutes until the solution becomes clear.
- 4. Add 200 μL Buffer F3 to the sample, mix completely by inverting the vial for 10 times. Incubate at room temperature for 2 min. Transfer the whole lysate to a lysate clearance column.

  Note: If processing more than 2 mL of culture, spin the lysate at 12,000 rpm for 2 min, transfer the relatively clear lysate to a lysate clearance column.
- Centrifuge at 12,000 rpm for 30 seconds.
   Note: If the lysate still remains in the column, spin for another 30 seconds.
- 6. Discard the lysate clearance column and add  $200 \,\mu\text{L}\ 100\%$  ethanol to the flow through in the collection tube, mix well by pipetting and transfer the sample to a DNA mini column.
- 7. Spin at 12,000 rpm for 30 seconds. Discard the flow through and reuse the collection tube.

- 8. Optional: Add 500 μL Buffer KB and centrifuge at 12,000 rpm for 30s. Discard the flow-through and put the column back to the collection tube.
  - Note: This step is NOT necessary if the plasmid is being purified from endA- strain such as DH5a and Top 10. Buffer KB wash is necessary for endA+ strains such as HB101, JM110, JM 101 and their derived strains.
- 9. Add 500 μL DNA Wash Buffer (Add ethanol to DNA wash buffer before use) into the spin column, centrifuge at 12,000 rpm (14,000 18,000 x g) for 30 seconds at RT. Remove the spin column from the tube and discard the flow-through.
- 10. Reinsert the spin column into the collection tube and centrifuge for 1 minute at 13,000 rpm. Note: Residual ethanol will be removed more effectively with the column lid open.
- 11. Carefully transfer the spin column into an Elution (collection) tube (Provided) and add 50-100 µL Elution Buffer into the column and let it stand for 1 minute. Elute the DNA by centrifugation at 12,000 rpm (14,000-18,000 x g) for 30 seconds to elute DNA.
  - Optional: The first elution normally yields around 70% of the plasmid DNA bound to the column. Add the eluted DNA back to the column for another elution yields 20-30% of the DNA bound to the column.