

Express Hi Fi DNA Polymerase (2 U/ μ L)

Introduction

To ensure successful PCR using **Express Hi Fi DNA Polymerase**, the following guidelines are provided to cover standard PCR. PCR of templates with high secondary structure, low template concentrations require further optimization.

Components

Cat#	EHF DNA Polymerase	5x Buffer	20 mM MgCl ₂	DMSO
EHF1101-02	2500 Units (1250 μ L)	10x1200 μ L	3 mL	1000 μ L
EHF1101-01	500 Units (250 μ L)	2x1200 μ L	1 mL	500 μ L

PCR Protocol

- Setup:** Assemble all reaction components on ice and quickly transferring the reactions to a thermocycler preheated to the denaturation temperature (95°C). All components should be mixed well and quickly spin down before start.

Component	20 μ l Reaction	50 μ l Reaction	Final Concentration
5X Express Hi Fi Buffer	4 μ l	10 μ l	1X
10 mM dNTPs	0.4 μ l	1 μ l	200 μ M
10 μ M fwd primer	1 μ l	2.5 μ l	0.5 μ M
10 μ M Rv Primer	1 μ l	2.5 μ l	0.5 μ M
Template DNA	variable	variable	< 250 ng
DMSO (optional)	(0.6 μ l)	(1.5 μ l)	3%
Express Hi Fi DNA Polymerase	0.5 μ l	1 μ l	1.0 units/50 μ l PCR
Nuclease-free water	to 20 μ l	to 50 μ l	

- Notes: Gently mix the reaction. Collect all liquid to the bottom of the tube by a quick spin if necessary. Overlay the sample with mineral oil if using a PCR machine without a heated lid.

Transfer PCR tubes from ice to a PCR machine with the block preheated to 98°C and begin thermocycling:

Thermocycling conditions for a routine PCR

STEP	TEMP	TIME
Initial Denaturation	95°C	30 seconds
25-35 Cycles	95°C 52-68°C 72°C	5-10 seconds 10-30 seconds 10 seconds per kb
Final Extension	72°C	8 minutes
Hold	4-8°C	

Important notes:

- Use of high-quality DNA templates greatly enhances the success of PCR. Recommended amounts of DNA template for a 50 µl reaction are as follows:

DNA	Amount
genomic	50 ng–250 ng
plasmid or viral	10 ng –50 ng

- The final concentration of each primer in a reaction using Express Hi Fi DNA Polymerase may be 0.2–1 µM.
- Mg⁺⁺ is critical to achieve optimal activity with Express Hi Fi DNA Polymerase. The final Mg⁺⁺ concentration in 1X Express Hi Fi HF Buffer is 1.5 mM. The optimal Mg⁺⁺ concentration is affected by dNTP concentration; the template being used and supplements that are added to the reaction. This can also be affected by the presence of EDTA. Mg⁺⁺ can be optimized in 0.2 mM increments using the MgCl₂ provided.
- Amplification of difficult template, such as GC-rich sequences or secondary structure, may be optimized by the presence of DMSO (included). A final concentration of 2% DMSO is recommended.
- The Express Hi Fi DNA Polymerase PCR products have blunt ends. Blunt-end cloning is recommended.