



The Inventor of EZgene™ Plasmid Purification System

96 Well RNA Isolation Kit

R6811

Handbook VER 2013.06

For in vitro research use only

If crystals form in buffers, warm up at 37°C to dissolve before use

Table of Contents

Introduction.....	2
Illustrated Protocol.....	3
Kit Contents/Storage and Stability.....	4
Preparing Reagents/Cleaning Plates.....	5
Total RNA purification Spin Protocol.....	5
Total RNA purification Vacuum/Spin Protocol.....	8
Optional DNase I Digestion Protocol.....	11
Troubleshooting Guide.....	12

Introduction

Genomic DNA contamination is a major challenge in the RNA purification process. The most common method to remove genomic DNA contaminants is to use DNase digestion.

The EZgene™ 96 HP RNA Isolation Kits are designed for fast isolation of total cellular RNA in high-throughput format without DNase I digestion. By using a special filter plate, genomic DNA can be effectively removed. The whole protocol can be completed in less than 40 minutes. RNA purified using the EZgene™ DNA/RNA method is ready for applications such as RT-PCR, qPCR, differential display, microarrays, etc.


Principle

The EZgene™ RNA Isolation kit combines the reversible binding properties of RNA technology with a specially designed buffer system which can effectively remove DNA before RNA isolation. Samples are first lysed and homogenized in a specially designed denaturing buffer (LCT), which immediately inhibits the activity of RNase and DNase. The lysate is then passed through a gDNA Clearance Plate which traps the genomic DNA. After adjusting the binding condition, the flow-through lysate that contains RNA are bound to the RNA plate. After three wash steps, purified RNA is eluted with RNase-free water.

Storage and Stability

All components of the 96 Total RNA Kit plus should be stored at 22°C-25°C. Under these conditions, RNA has successfully been purified and used for RT-PCR after 12 months of storage. Under cool ambient conditions, precipitation may form in Buffer LY and RB Wash Buffer. The crystals may be dissolved by heating the buffer at 37°C.

96 Well RNA Isolation Flowchart



Lyse Samples



Adjust binding condition
and load samples to RNA
plate



Wash 3x



Elute

Kit contents

Catalog#	R6818-00	R6818-01	R6818-02
96 RNA Plate	1	4	12
96 DNA Clearance Plate	1	4	12
2 mL Deep Well Plate	2	4	24
500 µL 96 Well Plate	1	4	12
Buffer LY	20 mL	125 mL	400 mL
RB Wash Buffer	70 mL	280 mL	850 mL
RNA Wash Buffer	40 mL	160 mL	2 x 200 mL
DEPC Water	10 mL	40 mL	120 mL
User Manual	1	1	1

96-Well Collection Plates (2 ml) are reusable; see Page 5 for instructions. Buffer LY and RB Wash Buffer contains a chaotropic salt. Use gloves and protective eye wear when handling this solution.

Preparing Reagents

Dilute **RNA Wash Buffer** with absolute ethanol (96-100%) as follows

Kit	Ethanol to be added
R6818-00	Add 160 mL absolute ethanol
R6818-01	Add 800 mL absolute ethanol
R6818-02	Add 800 mL absolute ethanol per bottle

Buffer LY: Add 20 µl β-mercaptoethanol per 1mL of Buffer LY before use.

Cleaning of 96-Well collection plates: If extra plates are needed, please call our customer service department for ordering information. To re-use the 96-Well collection plates, rinse them thoroughly with tap water, incubate overnight in 0.2M NaOH/1mM EDTA, rinse with distilled water and dry by air.

Spin Protocol

Materials to be provided by user

- 96%-100% ethanol
- 70% ethanol
- Multichannel pipette
- RNase-free filter pipette tips
- Reagent reservoirs for multichannel pipettes
- Centrifuge with rotor for 96-well plates
- Disposable latex gloves
- 2ml 96-well deep well plate
- Sealing film
- Swing Bucket Centrifuge capable of 5,000 x g and Adaptor for 96 Deep Well plates

Before Starting:

- Prepare **Buffer LY** and **RNA Wash Buffer** According to Preparing Reagents Section

Procedure:

1. A). LYSIS OF MONOLAYER CULTURED CELLS GROWN IN A MULTI-WELL TISSUE

Remove the medium in the culture plate by pipetting completely. Add **150 µL** of **Buffer LY** directly to each well. Mix thoroughly by pipetting up and down 10 times. Transfer the cell lysate into a new 2 ml deep-well plate (supplied).

B).LYSIS OF SUSPENSION CULTURED CELLS: Spin down up to 5×10^5 cells per well at 300 x g for 2 minutes at 4-20°C. Remove the medium completely by pipetting. Add **150 µL** of **Buffer LY** directly to each well. Mix thoroughly by pipetting up and down 10 times. Transfer the cell lysate into a new 2 ml deep-well plate (supplied).

Note: Add **20 µl β-mercaptoethanol** per **1 ml** of **Buffer LY** before use. The complete removal of supernatant is critical for the RNA isolation.

2. Seal the plate with a sealing film and shake vigorously for 1 min. Spin down briefly to avoid cross contamination.

3. Remove the sealing film and add 1 volume of (**150 µL**) of **70 % ethanol** to the sample. Seal the plate with the sealing film and mix thoroughly for 1 min by shaking. Spin down briefly to avoid cross contamination.

4. Place the 96 well RNA plate on top of a 2 ml deep-well plate and carefully transfer the entire sample from Step 3 (including any precipitate) to each well of the RNA plate.

Note: Pipet up and down for 5 times and then transfer to the sample/ethanol mix to the 96 well RNA plate.

5. Load the RNA plate /2 ml deep-well plate into a microplate holder, and place the whole assembly into the rotor bucket of the centrifuge. Spin at 5,000 x g for 5 minutes at room temperature. Discard the flow through.

6. Add **500 µl** of **RB Wash Buffer** directly into the each well of the RNA plate and centrifuge at 5,000 x g for 5 minutes at room temperature. Discard the flow through and reuse the 2 mL Deep Well Plate.

7. Add **600 µl** of **RNA Wash Buffer** to each well of the RNA plate. Centrifuge at 5,000 x g for 5 minutes at room temperature. Discard the flow-through and re-use the 2 ml deep-well Plate.

8. Add **600 µl** of **RNA Wash Buffer** to each well of the 96 RNA plate. Centrifuge at 5,000 x g for 10 -15 minutes at room temperature. The prolonged centrifugation is necessary to dry the RNA plate.

Note: It is very important to dry the RNA plate completely before the elution step to remove residual ethanol that might otherwise interfere with downstream applications.

9. Place the RNA plate on top of a 500 µl Elution Plate (supplied)

10. Add **75-100 µl** of **DEPC-treated water** to each well. Seal the plate with sealing film. Make sure to add water directly onto the center of RNA matrix. Incubate for 1 min at room temperature. Centrifuge at 5,000 x g for 5 minutes at room temperature to elute the RNA. Store RNA at -20°C.

Note: Elution volume can vary according to user preference.

Vacuum/Spin Protocol

Note that all centrifugation steps must be carried out at room temperature.

Materials to be provided by user

- 96%-100% ethanol
- Multichannel pipette
- RNase-free filter pipette tips
- Reagent reservoirs for multichannel pipettes
- Centrifuge with rotor for 96-well plates
- 2ml 96-well deep well plate
- Sealing films
- Swing Bucket Centrifuge capable of 5,000 x g and Adaptor for 96 Deep Well plates
- 96 well Vacuum Manifold

Before Starting:

Prepare Buffer LY and RNA Wash Buffer According to Preparing Reagents Section

Procedure:

1. A). LYSIS OF MONOLAYER CULTURED CELLS GROWN IN A MULTI-WELL TISSUE: Remove the medium in culture plate by pipetting. Add **150 µL** of **Buffer LY** directly to each well. Mix thoroughly by pipetting up and down 10 times.

B). LYSIS OF SUSPENSION CULTURED CELLS: Transfer aliquots of up to 5×10^5 cells into the wells of a 96-well microplate. Spin the plate at 300 x g for 5 minutes. Remove the medium completely by pipetting. Add **150 µL** of **Buffer LY** directly to each sample. Mix thoroughly by pipetting up and down 10 times.

Note: Add **20 µl β-mercaptoethanol** per **1 ml** of **Buffer LY** before use. The complete removal of supernatant is critical for the RNA isolation.

2. Add 1 volume of 100% ethanol to each well. Seal the plate with a sealing film and mix well by shaking for 1 min.

3. PREPARE THE VACUUM MANIFOLD: Place a 2 ml deep well plate or waste collection tray inside the vacuum manifold base. Place manifold's top section squarely over its base. Place the RNA plate on the manifold's top section, making sure the RNA plate is seated tightly on the rubber ring. Connect the vacuum manifold to the vacuum source. Keep the vacuum switch off.

4. Carefully transfer the entire sample from Step 2 to each well of RNA plate. Seal the plate with a sealing film and switch on the vacuum source. Apply vacuum until all the sample contents pass through the well membranes. Ventilate and turn off the vacuum.

Note: If some of the well is clogged, remove the plate and place it on top of a 2 ml deep-well plate. Centrifuge at 5000 x g for 5 minutes.

5. Remove the sealing film and add **500 µl** of **Buffer RB** to each well, seal the RNA plate with the sealing film and switch on the vacuum source. Apply vacuum until the entire liquid pass through the well membranes. Ventilate and turn off the vacuum.

6. Add **500 µL** of **RNA Wash Buffer** directly into each well of the RNA plate. Apply vacuum until all the liquid passes through the membranes. Switch off the vacuum, and ventilate the manifold. Repeat this step once.

7. Remove the RNA plate from the manifold and strike the bottom of the plate on a stack of paper towels. Repeat several times until there's no liquid released onto the paper towel.

8. Place the RNA plate on top of a 500 µL elution plate (Supplied) and add **75-100 µL DEPC-treated water** to each well of the RNA plate. Incubate at room temperature for 1 minute.

9. Spin at 5,000 x g for 5 min to elute RNA. Store RNA at -20°C.

9. Add **800 µl RNA Wash Buffer** to each well of the RNA plate and apply vacuum until all the liquid passes through the well membranes. Switch off the vacuum.

Note: RNA Wash Buffer must be diluted with absolute ethanol before use. Refer to label on bottle for directions.

10. Add **800 µl of RNA Wash Buffer** to each well of the of RNA plate and apply the vacuum until transfer is complete. Switch off the vacuum and ventilate the manifold.

11. Remove the RNA plate from the top plate of the vacuum manifold, and place the RNA plate on top of a 2 ml deep-well plate.

12. Centrifuge at 5000 x g for 5 minutes to dry the membrane.

13. Remove the RNA plate and place it on top of a new 500 µl collection plate (supplied with the kit).

14. Add **75-100 µL of DEPC-treated water** to each well. Seal the plate with a sealing film. Make sure to add water directly onto the center of RNA matrix. Incubate for 3 minutes at room temperature. Centrifuge at 5,000 x g for 5 minutes at room temperature to elute the RNA.

15. Reload the eluted RNA back to the RNA plate for a 2nd elution yields another 20-30% of the RNA. The first elution normally yields 60-70% of the RNA.

DNase I Digestion Protocol (Optional)

Since the DNA clearance plate eliminates most of the DNA, DNase I digestion is not necessary for most downstream applications. DNase I set could be bought from Biomiga.

Note: DNase I is very sensitive and prone to physical denaturing; do not vortex the DNase I mixture. Mix gently by inverting the tube. Prepare the fresh DNase I digestion mixture before beginning the RNA isolation procedure.

- Standard DNase buffers may not be compatible with Biomiga's DNase I Digestion Set.

1. Follow the standard protocol until the samples completely pass through the RNA Plate (Steps 1-6). Then complete the procedure using the following steps:

A) Add 300µl RB Wash Buffer to each well of the RNA Plate and centrifuge at 4,000 x g for 1 min.

B) For each RNA sample, prepare the DNase I digestion mixture as follows:

Buffer	Volume per Prep
DNase I Digestion Buffer	73.5 µL
RNase Free DNase I (20 Kunitz/µl)	1.5 µL
Total Volume	75 µL

C). Pipet 75 µL DNase I digestion mixture directly onto the surface of the membrane in each well of the RNA Plate. Be certain to pipet the mixture directly onto each membrane, as DNA digestion might not be complete if some of the mixture is retained on the walls or the O-rings of the RNA Plate.

D) Incubate at room temperature (15-30°C) for 15 minutes.

2. Continued proceed to Step 7 on page 6 (Spin version) or Step 9 on

Troubleshooting Guide

Please use this guide to troubleshoot any problems that may arise. For further assistance, please contact us at 858-597-0602.

Possible Problems and Suggestions

Problem	Cause	Solution
Little or no RNA	RNA remains on the column	Repeat elution step
		Membrane Pre-heat DEPC-water to 70°C
	Membrane is overloaded	Reduce quantity of starting material
Clogged Membrane	Incomplete homogenization	Completely homogenize sample
		Increase centrifugation time
		Reduce amount of starting material
Degraded RNA	Starting Tissue Problems	Freeze starting material quickly in liquid nitrogen
		Do not store tissue prior to extraction unless they are lysed first
	RNase contamination	Follow protocol closely, and work quickly
		Ensure not to introduce RNase during the procedure. Check buffers for RNase contamination
Problem in Downstream applications	Salt carry-over during elution	Ensure RNA Wash Buffer add ethanol as indicated on bottle.
		1 X RNA Wash Buffer must be stored and used at room temperature.
		Repeat wash with RNA Wash Buffer
DNA contamination	DNA contamination	Digest with RNase-free DNase I Solution
Low Abs ratios	RNA diluted in acidic buffer or water	DEPC-treated water is acidic and can dramatically lower Abs ₂₆₀ values. Use TE buffer to dilute RNA prior to spectrophotometric analysis.

Related Products

Catalog #	Product Name	Preps	Price \$
R6311-01	Tissue RNA kit	50	150.00
R6311-02	Tissue RNA kit	250	650.00
R1011-01	RNAsecure Solution	50 mL	45.00
R1011-02	RNAsecure Solution	250 mL	150.00
R6312-01	Tissue RNA midi kit	10	80.00
R6312-02	Tissue RNA midi kit	20	160.00
R6314-01	Tissue RNA maxi kit	10	120.00
R6314-02	Tissue RNA maxi kit	25	270.00
R6811-01	96-well tissue RNA kit	4x96	780.00
R6811-02	96-well tissue RNA kit	20x96	3300.00
R6411-01	Blood RNA mini kit	50	150.00
R6411-02	Blood RNA mini kit	250	680.00
R6812-01	96-well blood RNA kit	4x96	780.00
R6812-02	96-well blood RNA kit	20x96	3500.00

Limited Use and Warranty

This product is intended for in vitro research use only. Not for use in human.

This product is warranted to perform as described in its labeling and in Biomiga's literature when used in accordance with instructions. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Biomiga. Biomiga's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of Biomiga, to replace the products, Biomiga shall have no liability for any direct, indirect, consequential, or incidental damage arising out of the use, the results of use, or the inability to use it product.

For technology support or learn more product information, please visit our website at www.biomiga.com or contact us at (858) 597-0602.

Biomiga, Inc.

10637 Roselle Street, Suite C

San Diego, CA 92121

Tel: (858) 597-0602

Fax: (858) 538-1698

Email: info@biomiga.com

Visit our web at www.biomiga.com and learn more about Biomiga products.