



Catalog No. RT0201-02

Introduction

One-Step RT-PCR is designed for the reverse transcription (RT) and polymerase chain reaction (PCR) amplification of a specific target RNA from either total RNA or mRNA. This system uses a mixture of MMLV Reverse Transcriptase and Hot start *Taq* DNA polymerase in an optimized reaction buffer, and can detect RNA targets up to 4.0 kb. The amount of input total RNA can range from 0.01 pg to 2 µg. Sufficient reagents are provided for 100 amplification reactions of 50 µl each.

Kit Contents

Contents	Volume
Enzyme mix: MMLV RT/ Hot start <i>Taq</i> Mix	50 µl
2× Reaction mix (optimized buffer contains dNTP and Mg)	1.25 ml
5 mM Magnesium Sulfate	1.0 ml

Storage

Store at -20 °C, One-Step RT-PCR is stable for 2 year when stored at -20 °C.

Features

- Convenient one-tube setup
- Optimized RT-PCR buffer and robustness
- High sensitivity, specificity, and reproducibility

Recommended RT-PCR reaction assembly

The following protocol is suggested as a starting point.

Components	50 ul Rxn	Final Concentration
2× Reaction mix	25 ul	1×
Enzyme mix	1.0 ul	
Forward Primer (10 uM)	1.0 ul	200 nM
Reverse Primer (10 uM)	1.0 ul	200 nM
RNA template	× ul	Variable (0.01pg - 2ug)
Final Volume (ul)	50 ul	

* If needed, the magnesium concentration can further be optimized

1. Assemble the reaction on ice.
2. Program the thermal cycler so that cDNA synthesis is followed immediately by PCR amplification.
 - cDNA synthesis:** 1 cycle: 40–50 °C for 10 min
 - Denaturation:** 1 cycle: 94 °C for 2 min
 - PCR amplification:** 40 cycles:
 - 94 °C for 15 s
 - 55–65 °C for 30 s
 - 68 °C for 1 min/kb
3. Analyze RT-PCR amplified products by gel electrophoresis.