Yeast Plasmid Miniprep Kit

Introduction	2
Storage and Stability	2
Kit Contents	3
Before Starting	4
EZgene TM Yeast Plasmid Miniprep Spin protocol	5
EZgene TM Yeast Plasmid Miniprep Vacuum/Spin protocol	8
Trouble Shooting Guide	9
Related Products	11
Limited Use and Warranty	12

Introduction

EZgeneTM yeast plasmid kit is designed for rapid and reliable isolation of high-quality plasmid DNA from yeast cultures. Utilizing the reversible nucleic acid-binding properties of our matrix, the plasmid DNA is bound to the matrix while proteins and other unwanted impurities are eliminated by wash buffer. Pure DNA is then eluted. Purified DNA can be directly used in downstream applications such as PCR, restriction digestion, and Southern Blot.

The Yeast Plasmid Mini Kit combines the power of spin column technology with the lyticase, glass beads and alkaline-SDS lysis of yeast cells to yield high quality plasmid DNA in less than 1 hour. The mini spin columns facilitate the binding, washing, and elution steps, thus enabling multiple samples to be processed simultaneously. The actual plasmid yields depend on copy numbers, yeast strain, and conditions of growth. Because of low copy numbers, the maximum yield from 5 mL yeast culture is around 1 μ g.

This protocol has been successfully used to isolate autonomous plasmids from S. cerevisiae. As a modified alkaline lysis procedure, genomic DNA is virtually eliminated from the preparation. Note that all centrifugation steps should be carried out at room temperature.

Storage and Stability

All EZgeneTM Yeast Plasmid Mini Kit components are guaranteed for at least 12 months from the date of purchase when stored as follows: Buffer YP I/RNase A mix at 4°C, Lyticase at -20°C, all other components at room temperature.

Capacity

The binding capacity for the mini column is 40 μ g of plasmid DNA.

Kit Contents

Catalog#	YD1271-00	YD1271-01	YD1271-02
Preps	4	50	250
ezBind TM DNA Mini Columns	4	50	250
Buffer YPI	1.5 mL	15 mL	70 mL
Buffer YPII	1.5 mL	15 mL	70 mL
Buffer YPIII	2 mL	20 mL	100 mL
Buffer SE	3.0 mL	30 mL	135 mL
Buffer KB	3.0 mL	28 mL	135 mL
DNA Wash Buffer	2 mL	12 mL	50 mL
Glass beads	250 mg	2.7 g	13 g
Lyticase (units)	1100 U	11,000 U	55,000 U
RNase A	10 µL	50 µL	210 µL
Elution Buffer	600 µL	10 mL	30 mL
Manual	1	1	1

Materials to Be Provided by Users

- ☑ Tabletop micro-centrifuge and nuclease-free 1.5 mL tubes.
- \square Water bath set to 30 °C.
- \square Absolute ethanol (96%-100%).

Before Starting

Briefly examine this booklet and become familiar with each step. Prepare all components and have the necessary materials ready before starting.

Important

☑ Dilute **DNA Wash Buffer** with **absolute ethanol** as follows :

YD1271-00: Add 8 mL absolute ethanol

YD1271-01: Add 48 mL absolute ethanol

YD1271-02: Add 200 mL absolute ethanol per bottle

✓ Prepare a lyticase stock solution with Buffer SE and aliquot into adequate portions. Store each aliquot at -20 °C and thaw before use. Each sample will require 30 µL of this solution.

YD1271-00: Dissolve with 130 µL Buffer SE

YD1271-01: Dissolve with 1.6 mL Buffer SE

YD1271-02: Dissolve with 8 mL Buffer SE

 \square Add vial of RNase A to bottle of YPI and store at 4°C

EZgeneTM Yeast Plasmid Miniprep Spin Protocol

- 1. Inoculate **5 mL YDP medium** with yeast carrying desired plasmid and grow at 30 °C with agitation for 16-24 h.
- 2. Pellet 1-3 mL yeast culture (use $< 2 \ge 10^7$ cells) by centrifugation at 5,000 \times g for 5 min at room temperature.
- Discard medium and resuspend cells in 480 μL Buffer SE with 30 μL lyticase solution. Resuspend the pellet by vortexing at maxi speed for 1 min. Complete resuspension of cell pellet is vital of obtaining good yields. Incubate at 30 °C for at least 30 min.
- 4. Pellet spheroblasts by centrifuging at 4,000 x g for 5 min at room temperature. Discard the supernatant completely.
- 5. Resuspend the spheroblasts pellet with 250μ L Buffer YPI.
- 6. Add **50 mg glass beads** and vortex at max speed for 5 min. Let the sample stand to allow the beads to settle. Transfer the supernatant to a new 1.5 mL centrifuge tube.
- Add 250 μL Buffer YPII and mix by inverting and rotating the tube 4-6 times to obtain a cleared lysate. Incubate at room temperature for 5 min.

Avoid vigorous mixing as this will shear chromosomal DNA and lower plasmid purity. Store Buffer YPII tightly capped.

- 8. Add 350 μ L Buffer YP III and mix completely by sharp handshaking several times until a flocculent white precipitate forms. Centrifuge at 13,000 × g for 10 min at room temperature.
- 9. Carefully transfer the clear supernatant to a DNA mini column. Ensure that the pellet is not disturbed and that no cellular debris is carried over into the column. Centrifuge at $10,000 \times g$ for 30 seconds. Discard the flow-through and put the column back to the collection tube.
- 10. Add **300 \muL Buffer KB.** Centrifuge at **10,000** × **g** for 30 seconds. Discard the flow-through and put the column back to the collection tube.
- Add 700 μL DNA Wash Buffer. Centrifuge at 10,000 × g for 1 min. Discard flow-through.

DNA Wash Buffer is supplied as a concentrate and must be diluted with absolute ethanol according to the instructions one bottle or on Page 3.

- 12. **Optional**: Repeat step 11.
- 13. Centrifuge the empty column, with the lid open, for 2 min at **13,000 x g** to dry the matrix.

This step removes residual ethanol from the matrix.

14. Place column into a clean 1.5 mL microcentrifuge tube. Add 50-100 μL Elution Buffer(10 mM Tris-HCL, pH 8.5) to the column matrix, let it stand by for 1 min, centrifuge at 13,000 x g for 1 min to elute DNA.

Optional: Add the eluted DNA back to the column for a second elution yields 20-30% of additional DNA. The first elution normally yields 60-70% of the DNA.

15. Yield and quality of DNA: Determine the absorbance of an appropriate dilution (20-50 folds) of the sample at 260 nm and 280 nm. The DNA concentration is calculated as follows:

DNA concentration = $A_{260} \times 50 \times$ (**Dilution Factor**) µg / mL

EZgeneTM Yeast Plasmid Mini Vacuum/Spin Protocol

Carry out cell culture, lysis and neutralization as indicated in previous section (**Steps 1-8**). Instead of continuing with centrifugation, follow steps as below.

Please read through previous section of this manual before using this protocol.

- 1. Prepare the vacuum manifold according to manufacturer's instructions and connect the column to the manifold.
- 2. Load the clear supernatant from **Step 8** in page 5 to the column.
- 3. Switch on vacuum source to draw the sample through the column, and then turn off the vacuum.
- 4. Add 500 μ L Buffer KB to the column, draw the wash buffer through the column by turning on the vacuum source.
- 5. Wash the column by adding 700 μ L DNA Wash Buffer. Draw the wash buffer through the column by turning on the vacuum source.
- 6. (Optional) Repeat this step with another 700 μL DNA Wash Buffer.
- Assemble the column into a 2 mL collection tube and transfer the column to a microcentrifuge. Spin at maximal speed (13,000 x g) for 2 min to dry the column.
- Place the column in a clean 1.5 mL microcentrifuge tube and add 50-100 μL Elution Buffer (10 mM Tris-HCL, pH 8.5). Let the column stand by for 1 min at room temperature and centrifuge at 13,000 x g for 1 min to elute DNA.

Trouble Shooting Guide

Problems	Possible cause	Suggestions
Low DNA yield	Poor cell lysis	Do not use more than 5 mL (with high copy plasmids or 10 mL with low copy plasmids) culture with the basic protocol.
Low DNA yield	Cells may not be dispersed adequately	Completely disperse the cell suspension by vortexing after adding Buffer YPI. After adding Buffer YPII, mix completely to obtain a clear lysate.
Low DNA yield	Buffer YP II, if not tightly closed, may need to be replaced.	Prepare as follows: 0.2 N NaOH, 1% SDS.
Low DNA yield	Yeast culture overgrown or not fresh.	Do not incubate cultures for more than 24 hr at 30 °C. Storage of cultures for extended periods prior to plasmid isolation is detrimental.
Low DNA yield	Low copy number plasmid used	Increase culture volume to 10 mL and scale up buffer volume.
No DNA eluted	Extended centrifugation during elution step at higher than 13,000 x g. Matrix may be present in eluate and cause abnormal OD readings.	If the centrifugation speeds higher than specified, some matrix residues may be co-purified with the plasmid DNA, but it will not interfere with PCR or restriction digests. Centrifuge the samples at suggested speed.
NO DNA eluted	Incomplete mixing with Buffer YPI	Repeat the procedure, this time making sure to vortex the sample with Buffer YPIII immediately and completely.

	Insufficient mixing	Increase incubation time
	with Buffer YPII	with Buffer YPII. Ensure
		that no visible cell clumps
		remain.
	DNA Wash Buffer	Prepare DNA Wash Buffer
No DNA eluted	not diluted with	as instructed above.
	absolute ethanol.	
High molecular	Over mixing of cell	Do not vortex or mix
weight DNA	lysate upon	aggressively after adding
contamination of	addition of Buffer	Buffer YPII. Adequate
product	YPII.	mixing is obtained by
		simply inverting and
		rotating tube to cover walls
		with viscous lysate.
Optical densities	Trace contaminants	Make sure to wash column
do not agree with	eluted from column	as instructed.
DNA yield on agarose	increase A 260.	Alternatively, rely on
gel		agarose gel/ethidium
		bromide electrophoresis for quantization.
RNA visible on	RNase A not added	Add RNase A to Buffer
agarose gel	to Buffer YPI.	YPI.
Plasmid DNA	Ethanol not	Centrifuge column as
floats out of well while	completely	instructed.
loading	removed from	
agarose gel	column following	
_	wash steps.	

Related Products

Catalog #	Product Name	Preps	Price \$
PD1211-01	Plasmid mini kit	50	45.00
PD1211-02	Plasmid mini kit	250	220.00
PD1213-01	Plasmid mini kit II	50	60.00
PD1213-02	Plasmid mini kit II	250	250.00
PD1411-01	Plasmid midi kit	10	52.00
PD1411-02	Plasmid midi kit	25	125.00
PD1511-01	Plasmid maxi kit	10	120.00
PD1511-02	Plasmid maxi kit	25	250.00
PD1611-01	Plasmid ezFilter mega 3 kit	2	80.00
PD1611-02	Plasmid ezFilter mega 3 kit	10	380.00
PD1811-01	96-well plasmid mini kit	4x96	400.00
PD1811-02	96-well plasmid mini kit	20x96	2000.00
PD1311-01	BAC/PAC isolation kit	50	90.00
PD1311-02	BAC/PAC isolation kit	250	90.00
YD1281-01	96-well Yeast plasmid kit	4x96	400.00
YD1281-02	96-well Yeast plasmid kit	20x96	2000.00
YD1291-01	ezFilter 96-well Yeast	4x96	500.00
	plasmid kit		
YD1291-02	ezFilter 96-well Yeast	20x96	2300.00
	plasmid kit		

Limited Use and Warranty

This product is warranted to perform as described in its labeling and in Biomiga's literature when used in accordance with instructions. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Biomiga. Biomiga's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of Biomiga, to replace the products, Biomiga shall have no liability for any direct, indirect, consequential, or incidental damage arising out of the use, the results of use, or the inability to use it product.

For technical support or learn more product information, please contact us at (858) 603-3219 or visit our website at <u>www.biomiga.com</u>