




# UTI Panel 1, 2 and 3 Real Time PCR Detection



Channel	Bacteria
JOE	<i>Serratia marcescens</i>
FAM	<i>Klebsiella pneumoniae</i>
ROX	<i>Morganella morganii</i>
Cy5	<i>Mycoplasma hominis</i>



Channel	Bacteria
JOE	<i>C. glabrata</i>
FAM	<i>Enterobacter cloacae</i>
ROX	<i>Staphylococcus aureus</i>
Cy5	<i>Streptococcus agalactiae</i>



Channel	Bacteria
JOE	<i>Enterococcus faecalis</i>
Rox	<i>Proteus mirabilis</i>
Cy5	<i>C. albicans</i>
FAM	Human TERT(optional)