



UTI Panel 4 and 5 Real TIME PCR Detection



Channel	Bacteria
JOE	<i>pseudomonas aeruginosa</i>
FAM	<i>Ureaplasma parvum</i> serovar or <i>Ureaplasma urealyticum</i>
ROX	<i>Escherichia coli</i>
Cys	<i>Providencia stuartii</i>



Channel	Bacteria
JOE	<i>Staphylococcus saprophyticus</i>
FAM	<i>C. parapsilosis</i>
ROX	<i>Klebsiella oxytoca</i>
Cy5	<i>C. tropicalis</i>